

SHIFT LOCK SERIES

3160 SHIFT Clutch Lock



3160



3160P



- Smooth gear and bearing mechanism for infinite adjustment
- Can be used for laminated and thermoplastic sockets
- Easily compensate for the variance between pin and socket alignment when use with the 4 hole plate
- Quiet Operation
- 1 year warranty

P/N	Description	Weight Limit
3160	Shift Clutch Lock	136kg
3160P	Shift Clutch Lock w/Plate	136kg

3170 SHIFT Shuttle Lock



3170



3170P



- Serrated pin attachment
- Can be used for laminated and thermoplastic sockets
- Easily compensate for the variance between pin and socket alignment when used with the 4 hole plate
- 1 year warranty

P/N	Description	Weight Limit
3170	Shift Shuttle Lock	136kg
3170P	Shift Shuttle Lock w/Plate	136kg



The plate has three lock locations to allow for easy offset in test sockets or definitive sockets.



Room No. 206, Yangpyeong Bldg.,
37-2 Yangpyeongdong 2 ga,
Yeongdeungpo-gu,
SEOUL, 150-102, KOREA
Tel : 82-2-2636-8912~3
Fax : 82-2-2636-8914

USA:
418 E, Esther Place
Placentia, CA 92870
Tel: 1-714-524-0663

Europe:
Solamed GmbH
Bellmannskamp 4a,
21339 Lueneburg, Germany
Tel: 49-4131-266-544